


## Vacuum Experiment Plate 1003166 Vacuum Bell Jar 1020809

### Instruction sheet

09/16 ALF



### 1. Safety instructions

- Protect the surface of the vacuum experiment plate from moisture, chemicals and mechanical affects to guarantee that the system remains vacuum-tight.
- Defective vacuum bell jars could implode.
- Before conducting experiments check the vacuum bell jar for any damage (entrapped air does not endanger operating safety).

## 2. Description

Vacuum experiment plate and vacuum bell jar are used to set up a vacuum chamber for experiments in the coarse (low) and fine vacuum range.

### 2.1 Vacuum Experiment Plate

The vacuum experiment plate comprises a metal plate with rubber sealing disc mounted on a tripod and two cock valves used for hose connection on the pump side as well as for ventilation. A center bore hole with M12 thread is used to secure experiment equipment. A vacuum-tight power feed with 4-mm safety sockets and two cables with 4-mm safety plugs are available for the power supply.

### 2.2 Vacuum Bell Jar

Vacuum bell jar made of glass with grip knob and polished flange to be set on top of the vacuum experiment plate

## 3. Technical data

### Vacuum experiment plate

Diameter:	250 mm
Height:	90 mm
Current feed-through:	2-pole with 4-mm safety sockets
Power feed:	via 2 cables approx. 1 m in length with 4-mm safety plugs
Electrical limit specs.:	max. 48 V, max. 12 A
Vacuum connection:	2 hose nozzles 12 mm and 8 mm Ø

### Vacuum bell jar

Inner diameter::	190 mm
Height:	220 mm

## 4. Operation

To perform experiments the following equipment is also required:

1 Rotary-Vane Vacuum Pump, Two-Stage	1003317
1 Vacuum Hose 8 mm	1002619

- Before performing the experiment check for any damage to the vacuum bell jar.
- Make sure that the sealing disc as well as the polished edge of the bell jar are dust-free.
- At the start of the experiment, press the vacuum bell jar against the plate until air pressure provides sufficient force against the rubber seal, then let the bell jar go.
- If the handle of the valve points downwards, both cocks are closed.
- If the handle of the valve points in the direction of one cock, that cock is open and the opposite cock is closed.
- After performing the experiment ventilate the chamber with the pump switched off and the evacuation valve closed.